

MEMStepper_DualLockIn Tutorial

Version 1.0 Feb.19 2026
UC Berkeley-MIT Cao & Tang Lab




alexw25.us@berkeley.edu | alexw25.us@utexas.edu

UC Berkeley
Engineering



MEMS Test Platform

all-reconfigurable-quantum-materials-lab



UC Berkeley
Engineering


1-MEMS Basic Control

Direct Set Mode

Stepping Mode


Positioning (Scan) Mode

Positioning (Scan) Mode
with PID Control



1-MEMS Basic Control

1-1 Direct Set Mode



```
% Button pushed function: SetButton
function SetButtonPushed(app, event)
    if isa(app.LockIn, 'visa') && strcmp(app.LockIn.Status, 'open')
        fprintf(app.LockIn, sprintf('auxv 1,%f\n', app.V1Knob.Value/app.volt_amp_ratio));
        fprintf(app.LockIn, sprintf('auxv 2,%f\n', app.V2Knob.Value/app.volt_amp_ratio));
        fprintf(app.LockIn, sprintf('auxv 3,%f\n', app.V3Knob.Value/app.volt_amp_ratio));
        app.updateVoltages();
    end

    app.setSteppingEnable(false); % Stepping needs to be reinitialized
    app.setPositioningEnable(false); % Positioning needs to be reinitialized
end
```

1-MEMS Basic Control

1-2 Stepping Mode

% Button pushed function: InitButton

```
function InitButtonPushed(app, event)
    app.PhaseEditField.Value = 0;
    app.StepsEditField.Value = 0;
    app.setVoltages(app.DriveVoltageEditField.Value, app.PhaseEditField.Value);
    app.setSteppingEnable(true);
    if app.SyncCapFittingCheckBox.Value
        app.addCapFittingPoint();
    end
end
```

% Button pushed function: Minus18Button

```
function Minus18ButtonPushed(app, event)
    app.PhaseEditField.Value = mod(app.PhaseEditField.Value - 15, 360);
    app.StepsEditField.Value = app.StepsEditField.Value - 0.125;
    app.setVoltages(app.DriveVoltageEditField.Value, app.PhaseEditField.Value);
    if app.SyncCapFittingCheckBox.Value
        app.addCapFittingPoint();
    end
end
```

-1/8

```
function setVoltages(app, volt, phase)
    if ~isa(app.LockIn, 'visa') || ~strcmp(app.LockIn.Status, 'open')
        error('USB connection not established');
    end
    return;
    fprintf(app.LockIn, sprintf('auxv 1,%f\n', volt/app.volt_amp_ratio*MENStepper.vphase(phase/180*pi, 0)));
    fprintf(app.LockIn, sprintf('auxv 2,%f\n', volt/app.volt_amp_ratio*MENStepper.vphase(phase/180*pi, -2*pi/3)));
    fprintf(app.LockIn, sprintf('auxv 3,%f\n', volt/app.volt_amp_ratio*MENStepper.vphase(phase/180*pi, 2*pi/3)));
    app.updateVoltages();
end
```

1-3 Positioning (Scan) Mode



(2) Set Scan Params

(3) Cap Readout Channel

```
function GoTo(app, volt, target, microstep, speed)

    curr_step = app.StepsPhyEditField.Value;
    total_steps = round(abs(target-curr_step)*(microstep));
    wait_time = 1/(speed * microstep);
    step_inc = 1/microstep*sign(target-curr_step);

    for i=1:total_steps
        curr_step = curr_step + step_inc;
        phase = curr_step * 120;

        app.setVoltages(volt, phase);
        pause(wait_time);
        if isa(app.video, 'videoinput') && isrunning(app.video)
            trigger(app.video);
        end
        app.StepsPhyEditField.Value = curr_step;
        if app.StepsCapZero > -100
            app.StepsCapEditField.Value = curr_step - app.StepsCapZero;
        end
    end
end
```

2-MEMS PID Control

2-1 [Set Point] Positioning (Scan) Mode with PID Control

MEMS control & readout

SR830	SR844	Direct	Stepping	Positioning	Positioning (Feedback)	3.00	Wait Time (s)
GPIB0::8::...	GPIB0::8::...	Drive Voltage	50	Scan Step	0 ~ 7	Final Step	0
Connect	GPIB0::8::INSTR (SR830)	2	PID wait (s)	0.5	Log/Plot	Off <input type="checkbox"/> On <input type="checkbox"/>	Step Min.Loss
OUTPUT OFF	GPIB0::9::INSTR (SR844)	Init/Reset	Scan Microstep	2	PID Microstep	inf	Kp 0.35
							Ki 0.0005
							Kd 0.05
							Current Step
							Set Vd
							0.5
							Set Vz
							Zero Vz
							Flash Vz
							0.00

(1) Connect Lock-In First

LockIn % VISA object for SR830 (typically for driving/control)
 LockIn_Read % VISA object for secondary Lock-In (SR844) for sensing



```
>> MEMStepper_DualLockIn
Warning: No devices were detected for the 'winvideo' adaptor. For troubleshooting device detection issues, click here.
Warning: instrfind will be removed in a future release. For objects of type
- serialport use serialportfind
- tcpclient use tcpclientfind
- topserver use topserverfind
- udpport use udpportfind
- visadev use visadevfind
- aardvark use aardvarkfind
- ni845x use ni845xfind
- iodevice with LegacyMode=false use iodevicefind
Warning: visa will be removed in a future release. Use visadev instead.

Note: To use visadev with iodevice, create an iodevice object with name-value argument LegacyMode=false.
[LOG] SR830 successfully connected.
Resource: GPIB0::8::INSTR
Status: OPEN
Instrument ID: Stanford_Research_Systems,SR830,s/n22628,ver1.01
Warning: instrfind will be removed in a future release. For objects of type
- serialport use serialportfind
- tcpclient use tcpclientfind
- topserver use topserverfind
- udpport use udpportfind
- visadev use visadevfind
- aardvark use aardvarkfind
- ni845x use ni845xfind
- iodevice with LegacyMode=false use iodevicefind
Warning: visa will be removed in a future release. Use visadev instead.

Note: To use visadev with iodevice, create an iodevice object with name-value argument LegacyMode=false.
[LOG] SR844 successfully connected.
Resource: GPIB0::9::INSTR
Status: OPEN
Purpose: High-frequency Readout via OUTP? 1
Instrument ID: Stanford_Research_Systems,SR844,s/n43265,ver1.006
```

2-MEMS PID Control

2-1 [Set Point] Positioning (Scan) Mode with PID Control

(2) Set Params

- **Scan Step 0~7** Numeric Input start/stop step of calibration scan (physical step value)
- **Final Step** Input final target step (position after scan calibration)
- **Filter Window** Set the size of data filter window (number of data points). Used to smooth capacitance/step data (reduce noise) during **PID** positioning.
Recommended value: 1 (keep unchanged)
- **Speed (step/s)** Scan movement speed. unit: step/s, `wait_time = max(app.default_wait_time, 1/speed)`
- **PID wait (s)** PID sampling wait time input box. **Recommended value: 1/Speed**
- **Log/Plot Switch** Control real-time log plotting for PID feedback (capacitance-step curve). On: Plot curve + data points; Off: Only log data (no plot)
- **Step Min. Loss** This (step_tolerance) is the minimum step error threshold for PID control. It prevents integral saturation by stopping the accumulation of the **Integral (I)** term when the actual step error is smaller than this value. **Recommended value: 0.0001 to keep Integral (I) always "On"**
- **Scan Microstep** Select microstep resolution for scan (2/4/8...256; higher = more precise positioning). **Recommended value: 2 to keep linear**
- **PID Microstep** PID microstep resolution for feedback (2/4/8...256/**inf**; inf = continually moving). **Recommended value: inf**

2-MEMS PID Control

2-1 [Set Point] Positioning (Scan) Mode with PID Control

SR830

SR844

GPIB0::8:....

Connect

GPIB0::9:....

Connect

OUTPUT OFF

Direct

Stepping

Positioning

Positioning (Feedback)

Drive Voltage

50

Scan Step

0 ~ 7

Final Step

0

Filter Win.

1

Kp

0.35

CANCEL

0

Current Step

Speed (step/s)

2

PID wait (s)

0.5

Log/Plot

Off

On

Step Min.Loss

0.0001

Ki

0.0005

Kd

0.05

SCAN

GO

Init/Reset

Scan Microstep

2

PID Microstep

inf

0.00

Wait Time (s)

0.5

Set Vd


Set Vz

0.00

Zero Vz

Flash Vz

(3) Scan to Build Linear Model



```
[LOG] goTo_fit (scan version): wait_time=0.5000 s
[LOG] goTo_fit (scan version): wait_time=0.5000 s
[LOG] SCAN_PART Channel : R
[LOG] SCAN_PART_1 (Forward 15 points): k1 = 2.903e-05, b1 = 5.881e-04
[LOG] SCAN_PART_2 (Backward 15 points): k2 = 2.895e-05, b2 = 5.885e-04
[LOG] SCAN_TOTAL (Symmetric): k_total = 2.899e-05, b_total = 5.883e-04
[LOG] HYSTERESIS_AVG: Mean Delta Cap (dy) = 4.766e-07
[LOG] HYSTERESIS_AVG: Mean Delta Step (dx) = 0.0164 steps
(0.000000, 5.878450e-04)
(0.500000, 6.037570e-04)
(1.000000, 6.179570e-04)
(1.500000, 6.319560e-04)
(2.000000, 6.458550e-04)
(2.500000, 6.581410e-04)
(3.000000, 6.751610e-04)
(3.500000, 6.898650e-04)
(4.000000, 7.048700e-04)
(4.500000, 7.188690e-04)
(5.000000, 7.347810e-04)
(5.500000, 7.448520e-04)
(6.000000, 7.630810e-04)
(6.500000, 7.758710e-04)
(7.000000, 7.929920e-04)
(7.000000, 7.930920e-04)
(6.500000, 7.764750e-04)
(6.000000, 7.624770e-04)
(5.500000, 7.446510e-04)
(5.000000, 7.337740e-04)
(4.500000, 7.193730e-04)
(4.000000, 7.052730e-04)
(3.500000, 6.902670e-04)
(3.000000, 6.749590e-04)
(2.500000, 6.578390e-04)
(2.000000, 6.465590e-04)
(1.500000, 6.317550e-04)
(1.000000, 6.176550e-04)
(0.500000, 6.042610e-04)
(0.000000, 5.889530e-04)
```

Lower is better

2-MEMS PID Control

2-1 [Set Point] Positioning (Scan) Mode with PID Control

SR830

SR844

GPIB0::8:....

GPIB0::9:....

Connect

Connect

OUTPUT OFF

MEMS control & readout

Direct

Stepping

Positioning

Positioning (Feedback)

Drive Voltage

50

Scan Step

0 ~ 7

Final Step

0

Filter Win.

1

Kp

0.35

CANCEL

0

Current Step

Speed (step/s)

2

PID wait (s)

0.5

Log/Plot

Off

On

Step Min.Loss

0.0001

Ki

0.0005

Init/Reset

Scan Microstep

2

PID Microstep

inf

Kd

0.05

SCAN

GO

0.00

Wait Time (s)

0.5

Set Vd

Set Vz



Zero Vz

Flash Vz

0.00

(4) PID Positioning

```
Cap(theor-act): [ 5.919e-04 - 5.880e-04 = 3.891e-06]
Cap(theor-act): [ 5.919e-04 - 5.902e-04 = 1.776e-06]
Cap(theor-act): [ 5.919e-04 - 5.918e-04 = 1.639e-07]
Cap(theor-act): [ 5.919e-04 - 5.920e-04 = -3.709e-08]
Cap(theor-act): [ 5.919e-04 - 5.913e-04 = 6.679e-07]
Cap(theor-act): [ 5.919e-04 - 5.918e-04 = 1.639e-07]
Cap(theor-act): [ 5.919e-04 - 5.920e-04 = -3.709e-08]
Cap(theor-act): [ 5.919e-04 - 5.922e-04 = -2.381e-07]
Cap(theor-act): [ 5.919e-04 - 5.913e-04 = 6.679e-07]
Cap(theor-act): [ 5.919e-04 - 5.918e-04 = 1.639e-07]
Cap(theor-act): [ 5.919e-04 - 5.916e-04 = 3.659e-07]
Cap(theor-act): [ 5.919e-04 - 5.924e-04 = -4.401e-07]
Cap(theor-act): [ 5.919e-04 - 5.920e-04 = -3.709e-08]
Cap(theor-act): [ 5.919e-04 - 5.923e-04 = -3.391e-07]
Cap(theor-act): [ 5.919e-04 - 5.919e-04 = 6.391e-08]
Cap(theor-act): [ 5.919e-04 - 5.911e-04 = 8.699e-08]
Cap(theor-act): [ 5.919e-04 - 5.922e-04 = -2.381e-07]
Cap(theor-act): [ 5.919e-04 - 5.919e-04 = 6.391e-08]
Cap(theor-act): [ 5.919e-04 - 5.913e-04 = 6.679e-07]
Cap(theor-act): [ 5.919e-04 - 5.923e-04 = -3.391e-07]
Cap(theor-act): [ 5.919e-04 - 5.920e-04 = -3.709e-08]
Cap(theor-act): [ 5.919e-04 - 5.916e-04 = 3.659e-07]
>>
CWB(CPWB01-POS2) [ 2.9759e-04 - 2.9759e-04 = -3.6289e-06]
CWB(CPWB01-POS2) [ 2.9759e-04 - 2.9759e-04 = -3.7209e-06]
CWB(CPWB01-POS2) [ 2.9759e-04 - 2.9759e-04 = -3.7209e-06]
CWB(CPWB01-POS2) [ 2.9759e-04 - 2.9759e-04 = -3.7209e-06]
CWB(CPWB01-POS2) [ 2.9759e-04 - 2.9759e-04 = -3.7209e-06]
```




2-MEMS PID Control

2-3 [Scan Point] Positioning (Scan) Mode with PID Control

(1) Set Params

- **Speed (step/s)** defines the movement time of each step (eg. from “2” to “3”)
 - **PID wait time (s)** the delay for PID sampling
 - **Micro Step** Microstep count within each Step
- (2) Set switch**
- **Total time per step = 1/speed**
 - **Total time per micro step = 1/(speed*Micro Step)**
 - **Total PID sample count for each micro step = Total time per micro step / PID wait time**

(2) Set switch "On" to enable PID



MEMS control & readout

SR830 SR844

GPIO0::8::... GPIO0::9::...

Connect Connect

OUTPUT OFF

Direct Stepping Positioning Positioning (Feedback)

Drive Voltage 50 Target Step 2

Speed step/s 0.01

Micro Step 32

LogCapData Zero Cap

Steps (Phy) 0.2188

Set Vd 0.5

Set Vg

Init/Reset

GO

Close Loop Off On

(3) Start scanning: "GO"

(4) Set switch "Off" to terminate loop while PID running

(3) Start scanning: “GO”

(4) Set switch “Off” to terminate loop while PID running

2-MEMS PID Control


2-3 [Scan Point] Positioning (Scan) Mode with PID Control

Save “mems.log” for data plot/analysis

```
Cap(theor-act): [ 5.920e-04 - 5.920e-04 = 2.347e-09] | Step Loss (=Cap/k_fit): [ 0.000] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [ 0.0000 + -0.0000 + 0.0007
Cap(theor-act): [ 5.920e-04 - 5.921e-04 = -9.865e-08] | Step Loss (=Cap/k_fit): [ -0.003] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [-0.0013 + -0.0000 + -0.0007
Cap(theor-act): [ 5.920e-04 - 5.915e-04 = 5.063e-07] | Step Loss (=Cap/k_fit): [ 0.017] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [ 0.0066 + -0.0000 + 0.0041
Cap(theor-act): [ 5.920e-04 - 5.923e-04 = -2.997e-07] | Step Loss (=Cap/k_fit): [ -0.010] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [-0.0039 + -0.0000 + -0.0055
Cap(theor-act): [ 5.920e-04 - 5.921e-04 = -9.865e-08] | Step Loss (=Cap/k_fit): [ -0.003] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [-0.0013 + -0.0000 + 0.0014
Cap(theor-act): [ 5.920e-04 - 5.920e-04 = 2.347e-09] | Step Loss (=Cap/k_fit): [ 0.000] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [ 0.0000 + -0.0000 + 0.0007
Cap(theor-act): [ 5.920e-04 - 5.925e-04 = -5.017e-07] | Step Loss (=Cap/k_fit): [ -0.017] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [-0.0065 + -0.0000 + -0.0034
Cap(theor-act): [ 5.920e-04 - 5.913e-04 = 7.073e-07] | Step Loss (=Cap/k_fit): [ 0.024] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [ 0.0092 + -0.0000 + 0.0083
[LOG] Appended 200 points to capdata.log @ 2026-02-18 19:46:29.443288
[LOG] Logged successfully @ 2026-02-18 19:46:29.443288
[LOG] Feedback loop completed for step 7/64. Moving to next step.
```

```
[LOG] Starting feedback loop with total_cnt=200 for step 8/64 (scan_loop_step=2.500e-01).
```

```
Cap(theor-act): [ 5.929e-04 - 5.924e-04 = 5.150e-07] | Step Loss (=Cap/k_fit): [ 0.018] | Gains(P,I,D): [3.80e-01, 0.00e+00, 1.00e-01] | PID_Out(P+I+D): [ 0.0067 + -0.0000 + -0.0013
Cap(theor-act): [ 5.929e-04 - 5.947e-04 = -1.801e-06] |
Cap(theor-act): [ 5.929e-04 - 5.926e-04 = 3.140e-07] |
Cap(theor-act): [ 5.929e-04 - 5.939e-04 = -9.960e-07] |
Cap(theor-act): [ 5.929e-04 - 5.938e-04 = -8.950e-07] |
Cap(theor-act): [ 5.929e-04 - 5.926e-04 = 3.140e-07] |
Cap(theor-act): [ 5.929e-04 - 5.928e-04 = 1.120e-07] |
Cap(theor-act): [ 5.929e-04 - 5.925e-04 = 4.140e-07] |
Cap(theor-act): [ 5.929e-04 - 5.929e-04 = 1.202e-08] |
Cap(theor-act): [ 5.929e-04 - 5.933e-04 = -3.910e-07] |
Cap(theor-act): [ 5.929e-04 - 5.932e-04 = -2.910e-07] |
Cap(theor-act): [ 5.929e-04 - 5.918e-04 = 1.119e-06] |
Cap(theor-act): [ 5.929e-04 - 5.933e-04 = -3.910e-07] |
Cap(theor-act): [ 5.929e-04 - 5.932e-04 = -2.910e-07] |
```



mems.log

3-More

3-1

Script (Matlab R2024b)

[MEMStepper_DualLockIn.m](#)

App

(Matlab R2025b)

[MEMStepper_DualLockIn.mlapp](#)

```
% Create ImageMotionDetectionPanel
% app.ImageMotionDetectionPanel = uipanel(app.GridLayout);
% app.ImageMotionDetectionPanel.TitlePosition = 'centertop';
% app.ImageMotionDetectionPanel.Title = 'Image Motion Detection';
% app.ImageMotionDetectionPanel.Layout.Row = 2;
% app.ImageMotionDetectionPanel.Layout.Column = [3 4];
```

```
% [FIX] The modified code (the second section with [FIX] annotations) follows the constraint of MATLAB R2024 and ensures compatibility with both versions:
% 1. First, create the uipanel and directly bind it to the parent app.GridLayout at the time of, while setting basic properties simultaneously.
app.ImageMotionDetectionPanel = uipanel(app.GridLayout, 'Title', 'Image Motion Detection', 'BorderType', 'line', 'TitlePosition', 'centertop');
% 2. After the uipanel is successfully created and associated with its parent GridLayout
% set the Layout.Row and Layout.Column properties to define the component's position in the GridLayout.
app.ImageMotionDetectionPanel.Layout.Row = 2;
app.ImageMotionDetectionPanel.Layout.Column = [3 4];
```

3-2

Github

https://github.com/CHENfd515/MEMStepper_DualLockIn.git

3-2

Reference Article

[An adaptive moiré sensor for spectro-polarimetric hyperimaging](#)

